

EMERALD

SURGICAL INC.



ELECTROSURGICAL DEVICES

Emerald Surgical Inc.

528 Hewlett Street Franklin Square, NY 11010

www.emeraldsurgical.org

sales@emeraldsurgical.org

877-553-6372

ELECTROSURGICAL GENERATOR (OBS-350A)



Max Output Power : 350W

Intelligent Versatile Electrosurgical Generator
(Vessel Sealing) (Intelligent resistance feedback
automatic output compensation)

- ▶ Excellent compatibility
- ▶ Available in vessel sealing
- ▶ Suitable for general surgery, vascular surgery, hepatobiliary and pancreatic surgery, cardiothoracic surgery, gastroenterology, TURP and TUPV
- ▶ Compatible with argon coagulation, cystoscopy, arthroscopy, laparoscopy and hysteroscopy



Transport Data:

Packing box size	59*49.4*32(cm)
Gross weight	13.7 kg


Emerald Surgical, Inc.

528 Hewlett Street Franklin Square, NY 11010
www.emeraldsurgical.org - 877-553-6372
sales@emeraldsurgical.org

Single Use Sterile Disposable Hand Control Pencil (GBS-DB)

(Cautery Pencil)

Standard 3-pin connector with a 3mm cable
with blade electrode



Large color-coded
buttons with sealed
switching components



Ribbed electrode stem
provides easy insertion &
removal



Reusable Connecting Cable for Patient Plates (6070)

3 meter length for Valleylab



Blue Clip to connect
with Electrode



Valleylab REM



Emerald Surgical, Inc.


528 Hewlett Street Franklin Square, NY 11010
www.emeraldsurgical.org - 877-553-6372
sales@emeraldsurgical.org

FDA 510 (K)

Electrosurgical Grounding Pad

(GBS-DB 1031AC)

Electrosurgical Grounding Pad with Cable.




Foam backing	Adult	Split	Upright with Cable
The standard length of cable is 3 meter. Other length can be customized.			

Electrosurgical Grounding Pad

(GBS-DB 1031A)

Electrosurgical Grounding Pad without Cable.



Foam backing	Adult	Split	Upright
--------------	-------	-------	---------